




CLEARER CLUES FOR PERSISTENT DERMATITIS.

A ready-to-use core allergen patch test to assist qualified healthcare professionals in evaluating suspected allergic contact dermatitis.

2

READY-TO-USE
PANELS

19

ALLERGENS
AND MIXES


1

NEGATIVE
CONTROL

~48h

WEAR
PERIOD

READY TO APPLY. CLEARLY MAPPED.



PRE-COATED PANELS

Two individually
sealed panels, ready
for clinical use



PANEL-SPECIFIC MAPPING

Reading Cards align
reaction sites with
allergen positions



CLINICIAN-LED INTERPRETATION

Use with the current IFU,
documentation and clinical
judgement



APPLY



WEAR ~48 H



REMOVE



READ



CONTACT US ABOUT DERMAGY ACCESS, SAMPLES
AND LOCAL SUPPORT.



mag-pharma.com
info@mag-pharma.com

For professional use only. Refer to the current IFU for complete indications, contraindications, warnings and precautions.



DBS-20 BASELINE SERIES

19 ACTIVE ALLERGEN/MIX POSITIONS + 1 NEGATIVE CONTROL

PANEL A · POSITIONS 1–10

01	Nickel sulfate hexahydrate	2.5% pet.
02	Potassium dichromate	0.5% pet.
03	Cobaltous chloride hexahydrate	1.0% pet.
04	Fragrance mix I (incl. 8 allergens)	8.0% pet.
05	Fragrance mix II (incl. 6 allergens)	14.0% pet.
06	Peru balsam	25.0% pet.
07	Formaldehyde	2.0% aq.
08	Paraben mix (incl. 4 allergens)	16.0% pet.
09	Methylisothiazolinone + Methylchloroisothiazolinone	0.215% aq.
10	Quaternium-15	2.0% pet.

PANEL B · POSITIONS 11–20

11	Diazolidinyl urea	2.0% pet.
12	Thiuram mix (incl. 4 allergens)	1.0% pet.
13	Mercapto mix (incl. 4 allergens)	3.5% pet.
14	p-Phenylenediamine	1.0% pet.
15	Lanolin alcohols	30.0% pet.
16	Budesonide	0.01% pet.
17	Rosin	20.0% pet.
18	Oxybenzone	10.0% pet.
19	2-Hydroxyethyl methacrylate	2.0% pet.
20	Negative control	-

pet. = petrolatum (w/w) • aq. = aqueous/water (w/v)

EXTEND THE PANEL ACCORDING TO CLINICAL NEED



DME-20

MEDICAMENT & PRESERVATIVE

- Corticosteroids
- Topical antibiotics
- Antiseptics
- Preservatives



DEP-20

FRAGRANCE & PERSONAL CARE

- Fragrances
- Cosmetics & botanicals
- Cleansers & surfactants
- Product vehicles



DMO-20

MATERIALS & OCCUPATIONAL

- Resins & acrylates
- Rubber & PPE
- Hair chemicals
- Textile dyes

DISCUSS PANEL SELECTION AND AVAILABILITY WITH US.

mag-pharma.com info@mag-pharma.com

For professional use only. Refer to the current IFU for complete indications, contraindications, warnings and precautions.

